

ELECTRONIC CATALOGUE  
CUSTOMER'S AMENDMENTS  
ISSUE 1 — January 2012  
Page 1 of 16 pages

	<i>Amendment</i>																																												
	<p><b>AMENDMENTS</b></p> <table border="0"> <tr> <td>BA104-13 is now 58ml capacity</td> <td>BA105-11 is now 29ml capacity</td> </tr> <tr> <td>BA104-15 is now 118ml capacity</td> <td>BA105-13 is now 58ml capacity</td> </tr> <tr> <td>BA104-16 is now 207ml capacity</td> <td>BA105-15 is now 118ml capacity</td> </tr> <tr> <td>BA104-17 is now 532ml capacity</td> <td>BA105-17 is now 532ml capacity</td> </tr> <tr> <td>BA104-19 is now 710ml capacity</td> <td>BA105-19 is now 710ml capacity</td> </tr> <tr> <td>BA104-21 is now 710ml capacity</td> <td></td> </tr> <tr> <td>BA104-23 is now 798ml capacity</td> <td>BA107-15 is now 118ml capacity</td> </tr> <tr> <td>BA104-25 is now 1065ml capacity</td> <td>BA107-17 is now 532ml capacity</td> </tr> <tr> <td>BA104-27 is now 1242ml capacity</td> <td>BA107-19 is now 710ml capacity</td> </tr> </table> <p><b>NEW ITEMS</b></p> <p><b>Additional Whirl-Pak® bags.</b></p> <table border="0"> <tr> <td>BA105-16</td> <td>207ml, with writing panel, 180 x 95mm, pk.500</td> </tr> <tr> <td>BA105-31</td> <td>2721ml, with writing panel, 380 x 254mm, pk.250</td> </tr> <tr> <td>BA105-33</td> <td>3637ml, with writing panel, 508 x 254mm, pk.250</td> </tr> <tr> <td>BA106-15</td> <td>118ml, with writing panel and flat wire band, 185 x 75mm, pk.500</td> </tr> <tr> <td>BA107-25</td> <td>1065ml, with writing panel, self-standing, 380 x 125mm, pk.250</td> </tr> <tr> <td>BA107-27</td> <td>1242ml, with writing panel, self-standing, 380 x 150mm, pk.250</td> </tr> <tr> <td>BA107-29</td> <td>2041ml, with writing panel, self-standing, 380 x 190mm, pk.250</td> </tr> </table> <p><b><u>Wire racks, Whirl-Pak®, Nasco</u></b>  For convenient storage or transport of samples within Whirl-Pak bags. Dimensions shown are number of apertures x dimensions and to hold bags with capacities as indicated.</p> <table border="0"> <tr> <td>BA110-10</td> <td>21 x (40 x 50mm), up to 207ml bags</td> </tr> <tr> <td>BA110-20</td> <td>18 x (50 x 50mm), up to 207ml bags</td> </tr> <tr> <td>BA110-30</td> <td>20 x (64 x 64mm), 532 to 1065ml bags</td> </tr> <tr> <td>BA110-40</td> <td>6 x (114 x 114mm), 710 to 2041ml bags</td> </tr> <tr> <td>BA110-60</td> <td>Carrying handle for BA110-racks</td> </tr> </table> <p><b><u>Whirl-Pak® Scoop bags, Nasco</u></b>  Comprises 532ml Whirl-Pak bag, 230 x 115mm, containing a sterile, polystyrene, 60ml scoop. Supplied in a pack of 50.</p> <table border="0"> <tr> <td>BA115-10</td> <td>Scoop bags</td> </tr> </table>	BA104-13 is now 58ml capacity	BA105-11 is now 29ml capacity	BA104-15 is now 118ml capacity	BA105-13 is now 58ml capacity	BA104-16 is now 207ml capacity	BA105-15 is now 118ml capacity	BA104-17 is now 532ml capacity	BA105-17 is now 532ml capacity	BA104-19 is now 710ml capacity	BA105-19 is now 710ml capacity	BA104-21 is now 710ml capacity		BA104-23 is now 798ml capacity	BA107-15 is now 118ml capacity	BA104-25 is now 1065ml capacity	BA107-17 is now 532ml capacity	BA104-27 is now 1242ml capacity	BA107-19 is now 710ml capacity	BA105-16	207ml, with writing panel, 180 x 95mm, pk.500	BA105-31	2721ml, with writing panel, 380 x 254mm, pk.250	BA105-33	3637ml, with writing panel, 508 x 254mm, pk.250	BA106-15	118ml, with writing panel and flat wire band, 185 x 75mm, pk.500	BA107-25	1065ml, with writing panel, self-standing, 380 x 125mm, pk.250	BA107-27	1242ml, with writing panel, self-standing, 380 x 150mm, pk.250	BA107-29	2041ml, with writing panel, self-standing, 380 x 190mm, pk.250	BA110-10	21 x (40 x 50mm), up to 207ml bags	BA110-20	18 x (50 x 50mm), up to 207ml bags	BA110-30	20 x (64 x 64mm), 532 to 1065ml bags	BA110-40	6 x (114 x 114mm), 710 to 2041ml bags	BA110-60	Carrying handle for BA110-racks	BA115-10	Scoop bags
BA104-13 is now 58ml capacity	BA105-11 is now 29ml capacity																																												
BA104-15 is now 118ml capacity	BA105-13 is now 58ml capacity																																												
BA104-16 is now 207ml capacity	BA105-15 is now 118ml capacity																																												
BA104-17 is now 532ml capacity	BA105-17 is now 532ml capacity																																												
BA104-19 is now 710ml capacity	BA105-19 is now 710ml capacity																																												
BA104-21 is now 710ml capacity																																													
BA104-23 is now 798ml capacity	BA107-15 is now 118ml capacity																																												
BA104-25 is now 1065ml capacity	BA107-17 is now 532ml capacity																																												
BA104-27 is now 1242ml capacity	BA107-19 is now 710ml capacity																																												
BA105-16	207ml, with writing panel, 180 x 95mm, pk.500																																												
BA105-31	2721ml, with writing panel, 380 x 254mm, pk.250																																												
BA105-33	3637ml, with writing panel, 508 x 254mm, pk.250																																												
BA106-15	118ml, with writing panel and flat wire band, 185 x 75mm, pk.500																																												
BA107-25	1065ml, with writing panel, self-standing, 380 x 125mm, pk.250																																												
BA107-27	1242ml, with writing panel, self-standing, 380 x 150mm, pk.250																																												
BA107-29	2041ml, with writing panel, self-standing, 380 x 190mm, pk.250																																												
BA110-10	21 x (40 x 50mm), up to 207ml bags																																												
BA110-20	18 x (50 x 50mm), up to 207ml bags																																												
BA110-30	20 x (64 x 64mm), 532 to 1065ml bags																																												
BA110-40	6 x (114 x 114mm), 710 to 2041ml bags																																												
BA110-60	Carrying handle for BA110-racks																																												
BA115-10	Scoop bags																																												

**Biohazard Specimen bags, Nasco**

**Whirl-Pak<sup>®</sup> Trans-Pak<sup>®</sup> clinical specimen transport bags, non-sterile. Bag dimensions 230 x 150mm. With external document pouch 146 x 137mm and biohazard label. Accessory security tape is made of destructible material which cannot be removed in one piece and immediately tears if removal is attempted, indicating any tampering. Supplied in packs as indicated.**

**BA118-15 Biohazard specimen bags, pk.500**

**BA118-17 Security tape strips, pk.100**

**Soil Sample bags, Nasco**

**Used in sample submission for analytical assessment. Bleached paper with polyethylene liner and heat sealed base so that no sample sifting occurs. Capacity 769ml. With attached, folding tin strap tie. Supplied in a pack of 500.**

**BA120-25 Soil sample bags**

---

**CANCELLATIONS**

**BD160-02 Cancelled, offer BD160-05**

**BD166-10 Cancelled, offer BD166-11**

**NEW ITEM**

**SartoCollect data transfer software YSC02. Provides data transfer and integration into Microsoft<sup>®</sup> Excel, Word, Access or an ASCII file. Compatible with Windows 2000, XP, Vista and 7 operating systems.**

**BD166-11 YSC02**

---

**CANCELLATION**

**BH108-14 Cancelled, offer BH105-14**

**AMENDMENT**

**BH125-15 and -40 are now accessories for BH105 and BH120 series only**

**All other details remain unchanged.**

---

**CANCELLATION**

**BN240-08 Cancelled, no replacement**

**AMENDMENT**

**BS830-50 is PVC, not polyethylene as stated in the catalogue.**

**All other details remain unchanged.**

**CANCELLATIONS**

**CK324-15 Cancelled, no replacement**

**CK324-50 Cancelled, no replacement**

**CANCELLATION**

**CL867-07 Cancelled, offer CL867-06**

**AMENDMENT**

**CL867-06 is now for use with both CL865-20 and CL865-40 disinfectant gels.**

**All other details remain unchanged.**

**CANCELLATION**

**CX252-39 Cancelled, offer CX252-42**

**AMENDMENT**

**CX252-42 is now the lid for both CX250-31 and CX250-38**

**All other details remain unchanged.**

**CANCELLATION**

**DG100-10 Cancelled, offer DG100-05**

**NEW ITEM**

**Dialysis Tubing, Visking**

**Generally as DG100-10 but now supplied in a reel of 15 metres.**

**DG100-05 6 x 10mm, reel of 15m**

**CANCELLATION**

**DS605-10 Cancelled, offer DS605-15**

**NEW ITEM**

**Seeker**

**Generally as previous DS605-10 but with plastic handle.**

**DS605-15 Plastic**

**AMENDMENT**

**EL274-20 is no longer supplied with cooling gel pack.**

**All other details remain unchanged.**

**CANCELLATION**

**FC184-19 Cancelled, offer FC200-30**

**CANCELLATIONS**

**FZ530-10 Cancelled, no replacement**

**FZ530-20 Cancelled, no replacement**

**FZ530-30 Cancelled, no replacement**

**FZ530-40 Cancelled, no replacement**

**FZ530-50 Cancelled, no replacement**

**FZ530-60 Cancelled, no replacement**

**FZ533-10 Cancelled, no replacement**

**FZ533-30 Cancelled, no replacement**

**FZ534-10 Cancelled, no replacement**

**FZ534-20 Cancelled, no replacement**

**FZ534-30 Cancelled, no replacement**

**FZ534-40 Cancelled, no replacement**

**FZ534-50 Cancelled, no replacement**

**FZ534-60 Cancelled, no replacement**

**FZ537-10 Cancelled, no replacement**

**FZ537-20 Cancelled, no replacement**

**FZ537-30 Cancelled, no replacement**

**FZ537-40 Cancelled, no replacement**

**FZ537-50 Cancelled, no replacement**

**FZ537-60 Cancelled, no replacement**

**CANCELLATION**

**GM218-10 Cancelled, offer GM222-10**

**NEW ITEMS**

**Glass Cutting Knife**

**With plastic handle and replaceable, 40mm, tungsten-carbide blade.**

**GM222-10 Blade length 40mm**

**GM222-14 Spare blade**

---

**CANCELLATIONS**

**HE476-70 Cancelled, special order only**

**HE476-72 Cancelled, special order only**

---

**NEW ITEMS**

**Whirl-Pak<sup>®</sup> Sponge Probe, Nasco**

**Designed for surface sampling in difficult-to-reach areas. The sponge is adhered to a holder joint that is screwed onto a 160mm handle and ethylene oxide sterilised within a Whirl-Pak<sup>®</sup> 710ml bag with writing panel. In use a buffer is added to the bag and the excess is squeezed from the sponge which is then wiped over the test surface. The sponge is then replaced in the bag, the head is held through the bag walls and the handle unscrewed and discarded. The bag is then closed and the sponge sent for analysis. Supplied in a pack of 100.**

**HT850-15 Whirl-Pak<sup>®</sup> Sponge probes, pk.100**

**Poultry Rinse bags, Nasco**

**For testing E.coli presence in poultry meat. The sealed, sterile bag, 375 x 500mm W x D, has a perforation for easy opening and will accommodate small turkeys, chickens, and similar sized poultry. In use, the carcass is placed in the bag, buffer is added, and the bag is then twisted shut and folded around the bird. The bird is then rotated several times to rinse the carcass thoroughly, then the buffer liquid is removed for testing. Supplied in a pack of 250 bags.**

**HT880-10 Poultry rinse bags**

**NEW ITEMS**

**Slide mailers**

Generally as ML890-series, but with increased pack sizes.

**ML890-06** Single holder, pk.100

**ML890-11** Double holder, pk.100

**CANCELLATIONS**

**MN640-10** Cancelled, offer MN640-20

**MN640-25** Cancelled, offer MN640-20

**NEW ITEMS**

**Whirl-Pak<sup>®</sup> Blender Bags**

Addition to current MS865-range.

	Size, mm	Capacity, ml	Writing panel	Pack
<b>MS865-05</b>	150 x 230	710	No	<b>500</b>
<b>MS865-70</b>	185 x 330*	1538	Yes	<b>500</b>

\*With round bottom seal for use with circulating blenders.

**CANCELLATIONS**

**MW760-20** Cancelled, offer MW760-15

**MW760-30** Cancelled, offer MW760-25

## CANCELLATIONS

RG200-19 Cancelled, offer RG200-16  
RG200-42 Cancelled, offer RG200-44  
RG200-66 Cancelled, offer RG200-67  
RG200-87 Cancelled, offer RG200-85  
RG205-33 Cancelled, offer RG205-32

## NEW ITEMS

### Laboratory Refrigerators

Designed for the safe storage of chemicals, volatile liquids and specimens. Maintain internal air temperature between +2°C and +10°C.

- Automatic defrost
- Integral door lock and 2 keys
- Digital controller with LED temperature readout to 0.1oC
- Refrigerants:
  - R134a (RG200-16, -36 and -85 only)
  - R600a (RG200-44 and -67 only)
- “Spark-free” interior
- Adjustable shelf levels
- Audible and visual high and low temperature alarms
- Anti-microbial door handle

As described. With “spark-free” interior, automatic defrost and door lock. The thermostat is adjustable to maintain internal air temperature between +2°C and +10°C. All models are supplied with wire shelves and drawers as indicated. For 230V 50Hz single phase supplies.

	Model	Gross cap. litres	H x W x D mm	Shelves/drawers supplied
RG200-16	LR207C	82	660 x 502 x 540	2/1
RG200-36	LR307C	107	830 x 495 x 620	2/1
RG200-44	LSR151	151	845 X 558 X 600	3/1
RG200-67	LSR288	288	1570 x 595 x 634	5/1
RG200-85	LR1607C	475	1688 x 695 x 910	5/0

**Laboratory Refrigerators/Freezers**

Combined refrigerators/freezers for storage of temperature sensitive samples, volatile liquids, chemicals and poisons.

- Refrigerator range: +2°C to +10°C  
Freezer ranges: -6°C to -9°C (RG205-20)  
-18°C to -20°C (RG205-32, -47, -52)
- “Spark-free” interiors
- Integral door lock and 2 keys
- Digital controller with LED temperature readout
- Audible and visual high and low temperature alarms
- Adjustable self levels (refrigerator only)
- Anti-microbial door handle
- Refrigerants:
  - R134a (RG205-20 only)
  - R600a (RG205-32, -47 and -52)

**Laboratory Refrigerators/Freezers**

As described. With “spark-free” interior, door lock and wire shelves and drawers as indicated. For 230V 50Hz single phase supplies.

	Model	Gross cap.		H x W x D mm	Shelves/ drawers supplied
		refrigerator	freezer		
RG205-20	ISR27C	77	6.5	660 x 502 x 540	2/0
RG205-32	LSC119	106	13	860 x 540 x 600	2/1
RG205-47	LSC263	159	104	1680 x 540 x 600	3/2*
RG205-52	LSC324	162	162	1870 x 600 x 600	3/2*

\*RG205-47 and RG205-52 are supplied with half-width drawers and 3 and 4 freezer drawers, respectively.

**CANCELLATIONS**

- RG230-18 Cancelled, offer RG230-14
- RG230-23 Cancelled, offer RG230-22
- RG230-34 Cancelled, offer RG230-33
- RG230-44 Cancelled, offer RG230-43
- RG230-51 Cancelled, offer RG230-56
- RG230-59 Cancelled, offer RG230-58
- RG230-63 Cancelled, offer RG230-61
- RG230-72 Cancelled, offer RG230-71
- RG230-84 Cancelled, offer RG230-83
- RG230-96 Cancelled, offer RG230-94



**NEW ITEMS**

**Pharmacy Refrigerators**

Specifically designed for the storage of vaccines, insulin, temperature-sensitive samples and pharmaceuticals. Maintain internal air temperature between +2°C to +8°C.

- External digital temperature display with max./min. memory
- Forced air cooling for temperature stability and rapid temperature recovery after door opening
- Audio-visual alarm signal if temperature exceeds operating range. 107 litre and above models also incorporate a remote, auxiliary alarm output providing a mains failure alarm signal
- Integral door lock with 2 keys and open door alarm (except 41 and 45 litre models)
- Automatic defrost
- Refrigerant:                   - R134a (RG230-14, -16, -22, -26, -58, -61 and -64)  
  - R600a (RG230-33, -43, -56, -71, -83, -86 and -94)
- Adjustable shelf levels
- Built-in illumination and glass door models available
- Anti-microbial door handle

As described. With solid or glass door, door lock, adjustable wire shelves and wire base basket as indicated. For 230V 50Hz single phase supplies.

Model	Gross cap. litres	H x W x D mm	Shelves/ baskets supplied	Internal light
<b><u>Solid door models</u></b>				
RG230-14 PE102C	41	702 x 379 x 500	2/0	No
RG230-16 PE109C	45	530 x 470 x 465	1/1	No
RG230-22 PE207C	82	660 x 502 x 450	2/1	No
RG230-26 PE307C	107	830 x 495 x 620	2/1	No
RG230-33 PSR151	151	845 x 595 x 635	3/1	No
RG230-43 PSR273	273	1565 x 595 x 670	5/1	Yes
RG230-56 PSR334	334	1860 x 595 x 670	6/1	Yes
RG230-58 PE1607C	444	1688 x 685 x 910	5/0	No
<b><u>Glass door models</u></b>				
RG230-61 PG102C	41	702 x 379 x 500	2/0	No
RG230-64 PG207C	82	660 x 502 x 450	2/1	Yes
RG230-71 PGR151	151	845 x 595 x 635	3/1	Yes
RG230-83 PGR273	273	1565 x 595 x 670	5/1	Yes
RG230-86 PGR334	334	1860 x 595 x 670	6/1	Yes
RG230-94 PG1607C	444	1688 x 685 x 910	5/0	Yes

**NEW ITEM**

**Overboots**

Disposable, heavy duty polyethylene. With ankle ties. Will fit up to UK size 11 (EU size 46). Supplied in a pack of 50.

SA330-60 Overboots

---

**NEW ITEMS**

**Gloves, polyethylene**

SA350-25 Large, sterile, individually wrapped, pk.100 pairs

**Sleeves**

Disposable polyethylene, shoulder length gloves, size large. Supplied in a pack of 100.

SA360-10 Shoulder length gloves

---

**NEW ITEMS**

**Swing Sampler, Nasco**

Designed for sampling from a horizontal flowing stream. With 1.83m (6 foot) to 3.66m (12 foot) extending fibreglass pole with swing-hinged head and 960ml polyethylene bottle which is held firmly to the head with a snapper ring. The bottle is also supplied with a polyurethane lined cap for leakproof sample storage.

SD084-05 Swing sampler

**Tank Sampler, Nasco**

For liquid sampling at specific depths from tanks and vats. Comprises two 1.53m (5 foot) sections of 1.9cm ( $\frac{3}{4}$  inch) dia. plastic pipe that screw together and are marked at 0.3m (1 foot) increments. A valve at one end can be opened after lowering by pulling an internal stainless steel wire which is held whilst liquid flows in and released to shut the valve. The wire is pulled again or the bottom valve “tapped” to release the sample into a suitable container. The sampler holds approximately 89ml per 0.3m length and can be lengthened using an accessory 1.5m extension section.

SD084-10 Tank sampler, 3.05m

SD084-12 Accessory extension section, 1.5m

**Grease Trap Sampler, Nasco**

For measuring the grease level in a holding tank. Comprises two 1.22m (4 foot) transparent plastic tube sections that screw together and are 3.5cm in diameter to allow small food particles to pass easily through the valve opening at one end. The sections are marked at 0.3m (1 foot) intervals. The valve automatically opens when immersed and closes when the sample is removed. It can be re-opened by touching the valve top on a hard surface. Note: not suitable for solid or semi-solid materials.

SD084-33 Grease trap sampler

SD084-36 Cleaning brush for grease trap sampler

**Coliwasa Multi-phase Liquid Sampler, Nasco**

Enables representative sampling of multi-phase liquid wastes, operating within a wide range of viscous, corrosive, high volatiles or solids content liquids, except those containing ketones, nitrobenzene, dimethylformamide, mesityl oxide or tetrahydrofuran. Comprises a 1.53m (5 foot) x 1.9cm (<sup>3</sup>/<sub>4</sub> inch) i.d. tube, with valve at the base which is controlled by an internal rod. In use, the rod is pressed to close the valve, the tube is then immersed to the desired depth, then the valve opened by pulling the rod back 5 to 8cm. The rod is then pushed to close the valve and the sample removed.

SD084-44 Coliwasa multi-phase sampler

**Sampling Line, Nasco**

For use with Whirl-Pak<sup>®</sup> bags up to 2041ml capacity. Weighted urethane bag holder with strong clamp to secure the bag and 7.62m (25ft) nylon rope to enable sampling from liquid sources at deep levels. Supplied with carrying case and retaining clamps for different size bags. Weight 1.58kg.

SD084-55 Sampling line

**Fuel Sampler, Nasco**

For sampling ethanol and ethanol-blended fuels within petroleum tankers or fuel station forecourt tanks. Made up of three 1.53m (5 foot) clear PVC sections that screw together with a float valve at the base which opens on immersion and closes on withdrawal. Once the base of the tank has been reached, the sample can be withdrawn giving a core section of the fuel and any water layers present. The sampler is emptied by touching the valve tip against the tank side.

SD084-65 Fuel sampler

**NEW ITEMS**

**Additional accessories for Sludge Judge®**

**SD094-30** Cleaning rod, aluminium, to remove deposits.  
With slit to hold cleaning strips SD094-32

**SD094-32** Cotton cleaning strips for SD094-30, 50 x 300mm, pack of 50

**Sludge Judge II®, Nasco**

Generally as SD090-25 but with larger capacity. With 3.18cm (1¼ inch) diameter tubing, holding approximately 178ml (6oz) of liquid between the 0.30m (1ft) marks.

**SD095-15** Sludge Judge II®

**SD096-12** Accessory centre section, 1.53m (5ft)

**SD096-17** Spare top section with rope, 1.53m (5ft)

**SD096-22** Spare bottom section with valve

**Sludge Judge® Ultra, Nasco**

Generally as SD090-25 but constructed in a strong polycarbonate treated with UV stabiliser to reduce deterioration from exposure to the sun. Very rigid, with minimal flexing when filled and temperature resistant up to 138°C.

**SD097-20** Sludge Judge® Ultra

**SD098-16** Accessory centre section, 1.53m (5ft)

**SD098-18** Spare top section with rope, 1.53m (5ft)

**SD098-21** Spare bottom section with valve

**Dry Depth Sampler, Nasco**

Ideal for collecting dry, granular, free-flowing materials such as grain, corn or small pellets, at a specific depth. The sample head is screwed onto a 1.83m (6ft) to 3.05m (10ft) extendible aluminium pole. In use, the sampler is pushed down into the material to the depth required and the material flows through the head. When the head is withdrawn, an end flap shuts and the sample at that depth is retained.

**SD099-10** Dry depth sampler

---

**NEW ITEM**

**Scoops, polystyrene, sterile**

60ml capacity. One-piece moulded, with no seams or folds to retain sample. 200mm long. Individually wrapped within a pack of 100.

**SH207-30** Scoops, sterile, 60ml pk.100

**NEW ITEMS**

**Spatulas, sterile, disposable**

Polystyrene. Spoon one end, 14ml capacity. Individually bagged in an outer pack of 200.

SH420-10 178mm long, pk.200

**Spoons and Sporks, Nasco**

Disposable, sterile, white polystyrene. Supplied in sleeves of 20 within a pack of 100.

SH425-25 Spoons, pk.100

SH425-40 Sporks, pk.100

---

**CANCELLATIONS**

TP466-22 Cancelled, offer RG416-90

TP466-38 Cancelled, no replacement

TP466-42 Cancelled, offer RG416-95

---

**NEW ITEM**

**Settleometer Kit, Nasco**

For the evaluation of sludge quality and settlement. Comprises a 1.4 litre acrylic measuring cylinder, dual scale graduated in % and settled sludge volume (SSV) and supplied with a mixing paddle, snap-on lid and instructions. Overall H x dia. 216 x 114mm. Weight 0.6kg.

WR170-05 Settleometer kit

**NEW ITEMS**

**Thio-Bags<sup>®</sup>, Nasco**

**Generally as WR290-series but non self-standing.**

**WR290-02 Thio-Bags<sup>®</sup>, 100ml, pk.100**

**WR290-11 Thio-Bags<sup>®</sup>, 300ml, pk.100**

**WR290-95 Mailing tube for 100ml Thio-Bags<sup>®</sup>, rigid polyethylene,  
for multiple re-use**

**Racks for Thio-Bags<sup>®</sup>**

**Dimensions shown are number of apertures x dimensions of apertures.**

**BA110-10 21 x (40 x 50mm), for WR290-02 and -07**

**BA110-20 18 x (50 x 50mm), for WR290-02 and -07**

**BA110-30 20 x (64 x 64mm), for WR290-11 and -12**

***T H E E N D***

**PRICES FOR NEW ITEMS**

Page 15 of 16 Pages

January 2012

CAT. NO.	PRICE £	Per	CODE
BA105-16	54.50	pack	AW
BA105-31	120.00	pack	AW
BA105-33	140.00	pack	AW
BA106-15	39.65	pack	AW
BA107-25	80.00	pack	AW
BA107-27	90.00	pack	AW
BA107-29	110.00	pack	AW
BA110-10	23.00	each	AW
BA110-20	30.00	each	AW
BA110-30	30.00	each	AW
BA110-40	33.00	each	AW
BA110-60	17.00	each	AW
BA115-10	79.90	pack	AW
BA118-15	93.00	pack	AW
BA118-17	12.50	pack	AW
BA120-25	119.00	pack	AW
BD166-11	436.00	each	FO
DG100-05	89.25	each	AO
DS605-15	1.00	each	RO
	0.43	each per 10	
GM222-10	39.00	each	AO
GM222-14	22.30	each	AO
HT850-15	99.00	pack	AW
HT880-10	52.00	pack	AW
ML890-06	29.90	pack	AO
ML890-11	32.20	pack	AO
MS865-05	150.00	pack	AW
MS865-70	220.00	pack	AW
RG200-16	468.00	each	PO
RG200-36	540.00	each	PO
RG200-44	475.00	each	PO
RG200-67	675.00	each	PO
RG200-85	2848.00	each	PO
RG205-20	468.00	each	PO
RG205-32	375.00	each	PO
RG205-47	725.00	each	PO
RG205-52	1050.00	each	PO
RG230-14	585.00	each	PO
RG230-16	441.00	each	PO
RG230-22	754.00	each	PO

<b>RG230-26</b>	<b>980.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-33</b>	<b>895.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-43</b>	<b>1195.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-56</b>	<b>1795.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-58</b>	<b>2975.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-61</b>	<b>695.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-64</b>	<b>999.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-71</b>	<b>1195.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-83</b>	<b>1595.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-86</b>	<b>2195.00</b>	<b>each</b>	<b>PO</b>
<b>RG230-94</b>	<b>3229.00</b>	<b>each</b>	<b>PO</b>
<b>SA330-60</b>	<b>12.50</b>	<b>pack</b>	<b>AW</b>
<b>SA350-25</b>	<b>75.00</b>	<b>pack</b>	<b>AW</b>
<b>SA360-10</b>	<b>10.00</b>	<b>pack</b>	<b>AW</b>
<b>SD084-05</b>	<b>99.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-10</b>	<b>80.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-12</b>	<b>26.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-33</b>	<b>120.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-36</b>	<b>14.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-44</b>	<b>40.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-55</b>	<b>85.00</b>	<b>each</b>	<b>EO</b>
<b>SD084-65</b>	<b>90.00</b>	<b>each</b>	<b>EO</b>
<b>SD094-30</b>	<b>25.00</b>	<b>each</b>	<b>EO</b>
<b>SD094-32</b>	<b>10.00</b>	<b>pack</b>	<b>EO</b>
<b>SD095-15</b>	<b>110.00</b>	<b>each</b>	<b>EO</b>
<b>SD096-12</b>	<b>40.00</b>	<b>each</b>	<b>EO</b>
<b>SD096-17</b>	<b>28.00</b>	<b>each</b>	<b>EO</b>
<b>SD096-22</b>	<b>45.00</b>	<b>each</b>	<b>EO</b>
<b>SD097-20</b>	<b>99.00</b>	<b>each</b>	<b>EO</b>
<b>SD098-16</b>	<b>33.00</b>	<b>each</b>	<b>EO</b>
<b>SD098-18</b>	<b>34.00</b>	<b>each</b>	<b>EO</b>
<b>SD098-21</b>	<b>47.00</b>	<b>each</b>	<b>EO</b>
<b>SD099-10</b>	<b>95.00</b>	<b>each</b>	<b>EO</b>
<b>SH207-30</b>	<b>115.00</b>	<b>pack</b>	<b>AW</b>
<b>SH420-10</b>	<b>175.00</b>	<b>pack</b>	<b>AW</b>
<b>SH425-25</b>	<b>12.00</b>	<b>pack</b>	<b>AW</b>
<b>SH425-40</b>	<b>9.00</b>	<b>pack</b>	<b>AW</b>
<b>WR170-05</b>	<b>60.00</b>	<b>each</b>	<b>AW</b>
<b>WR290-02</b>	<b>17.50</b>	<b>pack</b>	<b>AW</b>
<b>WR290-11</b>	<b>71.00</b>	<b>pack</b>	<b>AW</b>
<b>WR290-95</b>	<b>4.00</b>	<b>each</b>	<b>AW</b>

***T H E E N D***